

anti- SPRY2 antibody

Product Information

Catalog No.:	FNab08199
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May function as an antagonist of fibroblast growth factor(FGF) pathways and may negatively modulate respiratory organogenesis.

Immunogen information

Immunogen:	sprouty homolog 2(Drosophila)
Synonyms:	hSPRY2, Protein sprouty homolog 2, sprouty homolog 2(Drosophila), Spry 2, SPRY2
Observed MW:	39 kDa, 65 kDa
Uniprot ID :	O43597

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:200-1:2000 IHC: 1:20-1:200 IF: 1:10-1:100
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using FNab08199 (SPRY2 Antibody) at dilution of 1:50



A375 cells were subjected to SDS PAGE followed by western blot with FNab08199 (SPRY2 antibody) at dilution of 1:500