

TRIM51 antibody

Product Information

Catalog No.:	FNab08202
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	SPRY domain containing 5
Synonyms:	Tripartite motif-containing protein 51 SPRY domain-containing protein 5 TRIM51 SPRYD5
Observed MW:	35 kDa
Uniprot ID :	Q9BSJ1

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab08202(SPRYD5 Antibody) at dilution of 1:300