

SDCBP2 antibody

Product Information

Catalog No.:	FNab08457
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	syndecan binding protein(syntenin) 2
Synonyms:	Syntenin-2 Similar to TACIP18 (SITAC) Syndecan-binding protein 2 SDCBP2 SITAC18
Observed MW:	35-37 kDa
Uniprot ID :	Q9H190

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab08457(SDCBP2 Antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab08457(Syntenin 2 antibody) at dilution of 1:1000