

TBX22 antibody

Product Information

Catalog No.:	FNab08539
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probable transcriptional regulator involved in developmental processes. This is major determinant crucial to palatogenesis.

Immunogen information

Immunogen:	T-box 22
Synonyms:	T-box transcription factor TBX22 (T-box protein 22) TBX22 TBOX22
Observed MW:	45 kDa and 55 kDa
Uniprot ID :	Q9Y458

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab08539(TBX22 Antibody) at dilution of 1:1000