

TCTN2 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab08568 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


Component of the tectonic-like complex, a complex localized at the transition zone of primary cilia and acting as a barrier that prevents diffusion of transmembrane proteins between the cilia and plasma membranes. Required for hedgehog signaling transduction(By similarity).

Immunogen information


| | |
|--------------|---------------------------------|
| Immunogen: | tectonic family member 2 |
| Synonyms: | Tectonic-2 TCTN2 C12orf38 TECT2 |
| Observed MW: | 77 kDa |
| Uniprot ID : | Q96GX1 |

Application

| | |
|-----------------------|---|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, IP, IF, IHC |
| Recommended dilution: | IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200 |
| Image: | |



Immunohistochemical of paraffin-embedded human liver using FNa08568(TCTN2 antibody) at dilution of 1:100



IP Result of anti-TCTN2 (IP: FNa08568, 4ug; Detection: FNa08568 1:500) with mouse lung tissue lysate 4000ug.