

## TFAP4 antibody

### Product Information

Catalog No.:	FNab08621
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Transcription factor that activates both viral and cellular genes by binding to the symmetrical DNA sequence 5'-CAGCTG-3'.

### Immunogen information


Immunogen:	transcription factor AP-4(activating enhancer binding protein 4)
Synonyms:	Transcription factor AP-4 Activating enhancer-binding protein 4 Class C basic helix-loop-helix protein 41 (bHLHc41) TFAP4 BHLHC41
Observed MW:	39 kDa
Uniprot ID :	Q01664

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab08621(TFAP4 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab08621(TFAP4 antibody) at dilution of 1:500