

TMEM176A antibody

Product Information

Catalog No.:	FNab08761
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	transmembrane protein 176A
Synonyms:	Transmembrane protein 176A Hepatocellular carcinoma-associated antigen 112 TMEM176A HCA112
Observed MW:	28 kDa
Uniprot ID :	Q96HP8

Application

Reactivity:	Human
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab08761(TMEM176A antibody) at dilution of 1:500