

C2CD2L antibody

Product Information

Catalog No.:	FNab08774
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	C2CD2-like
Synonyms:	Phospholipid transfer protein C2CD2L C2 domain-containing protein 2-like (C2CD2-like) Transmembrane protein 24 C2CD2L KIAA0285 TMEM24
Observed MW:	90 kDa
Uniprot ID :	O14523

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08774(TMEM24 antibody) at dilution of 1:1000