

## CLTRN antibody

### Product Information

Catalog No.:	FNab08778
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Regulator of SNARE complex function. Stimulator of beta cell replication.

### Immunogen information

Immunogen:	transmembrane protein 27
Synonyms:	Collectrin Transmembrane protein 27 CLTRN TMEM27
Observed MW:	32 kDa
Uniprot ID :	Q9HBJ8

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab08778(TMEM27 antibody) at dilution of 1:500