

TMEM70 antibody

Product Information

Catalog No.:	FNab08789
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4F2
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in biogenesis of mitochondrial ATP synthase.

Immunogen information


Immunogen:	transmembrane protein 70
Synonyms:	Transmembrane protein 70, mitochondrial TMEM70
Observed MW:	18 kDa
Uniprot ID :	Q9BUB7

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab08789(TMEM70 antibody) at dilution of 1:50



A2780 cells were subjected to SDS PAGE followed by western blot with FNab08789(TMEM70 antibody) at dilution of 1:500