

TPH1 antibody

Product Information

Catalog No.:	FNab08884
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	tryptophan hydroxylase 1
Synonyms:	Tryptophan 5-hydroxylase 1 Tryptophan 5-monooxygenase 1 TPH1 TPH TPRH TRPH
Observed MW:	50 kDa
Uniprot ID :	P17752

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HL60 cells were subjected to SDS PAGE followed by western blot with FNab08884(TPH1 antibody) at dilution of 1:1000