

## TPPP2 antibody

### Product Information

Catalog No.:	FNab08899
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May bind tubulin but has no microtubule bundling activity.

### Immunogen information

Immunogen:	tubulin polymerization-promoting protein family member 2
Synonyms:	Tubulin polymerization-promoting protein family member 2 Protein p25-beta TPPP/p18 TPPP2 C14orf8
Observed MW:	18 kDa
Uniprot ID :	P59282

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab08899(TPPP2 antibody) at dilution of 1:800