

TKT antibody

Product Information

Catalog No.:	FNab08938
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


This gene encodes a thiamine-dependent enzyme which plays a role in the channeling of excess sugar phosphates to glycolysis in the pentose phosphate pathway. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Immunogen information


Immunogen:	transketolase
Synonyms:	Transketolase (TK) TKT
Observed MW:	68 kDa
Uniprot ID :	P29401

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded rat brain using FNab08938(TKT antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab08938(TKT antibody) at dilution of 1:1000