

TSR2 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab09071 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


May be involved in 20S pre-rRNA processing.

Immunogen information


| | |
|--------------|---|
| Immunogen: | TSR2, 20S rRNA accumulation, homolog |
| Synonyms: | Pre-rRNA-processing protein TSR2 homolog TSR2 |
| Observed MW: | 21 kDa |
| Uniprot ID : | Q969E8 |

Application

| | |
|-----------------------|--|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IHC, IF |
| Recommended dilution: | WB: 1:500-1:5000; IHC: 1:20-1:200;IF: 1:10-1:100 |
| Image: | |



Immunohistochemistry of paraffin-embedded human liver using FNa09071(TSR2 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNa09071(TSR2 antibody) at dilution of 1:400