

TTC39B antibody

Product Information

Catalog No.:	FNab09086
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	tetratricopeptide repeat domain 39B
Synonyms:	Tetratricopeptide repeat protein 39B (TPR repeat protein 39B) TTC39B C9orf52
Observed MW:	66 kDa, 59 kDa
Uniprot ID :	Q5VTQ0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab09086(TTC39B antibody) at dilution of 1:800