

TXNL4B antibody

Product Information

Catalog No.:	FNab09131
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Essential role in pre-mRNA splicing. Required in cell cycle progression for S/G(2) transition.

Immunogen information

Immunogen:	thioredoxin-like 4B
Synonyms:	Thioredoxin-like protein 4B Dim1-like protein TXNL4B DIM2 DLP
Observed MW:	17 kDa
Uniprot ID :	Q9NX01

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09131(TXNL4B antibody) at dilution of 1:500