

UAP1L1 antibody

Product Information

Catalog No.:	FNab09141
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	UDP-N-acetylglucosamine pyrophosphorylase 1-like 1
Synonyms:	UDP-N-acetylhexosamine pyrophosphorylase-like protein 1 UAP1L1
Observed MW:	57 kDa
Uniprot ID :	Q3KQV9

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:20-1:200
Image:	



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with FNab09141(UAP1L1 Antibody) at dilution of 1:600