

## UBXN10 antibody

### Product Information

Catalog No.:	FNab09215
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	UBX domain protein 10
Synonyms:	UBX domain-containing protein 10 UBX domain-containing protein 3 UBXN10 UBXD3
Observed MW:	30-40 kDa
Uniprot ID :	Q96LJ8

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab09215(UBXN10 Antibody) at dilution of 1:1000