

## HAUS7 antibody

### Product Information

Catalog No.:	FNab09221
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	UCHL5 interacting protein
Synonyms:	HAUS augmin-like complex subunit 7 26S proteasome-associated UCH37-interacting protein 1 UCHL5-interacting protein X-linked protein STS1769 HAUS7 UCHL5IP UIP1
Observed MW:	46 kDa
Uniprot ID :	Q99871

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:1000; IHC:1:20-1:200;IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab09221(UCHL5IP Antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab09221(UCHL5IP Antibody) at dilution of 1:300