

## UTP23 antibody

### Product Information

Catalog No.:	FNab09349
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in rRNA-processing and ribosome biogenesis.

### Immunogen information

Immunogen:	UTP23, small subunit(SSU) processome component, homolog(yeast)
Synonyms:	rRNA-processing protein UTP23 homolog UTP23 C8orf53
Observed MW:	17 kDa
Uniprot ID :	Q9BRU9

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:5000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab09349(UTP23 antibody) at dilution of 1:1000