

## DCAF4L2 antibody

### Product Information

Catalog No.:	FNab09487
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	WD repeat domain 21C
Synonyms:	DDB1- and CUL4-associated factor 4-like protein 2 WD repeat-containing protein 21C DCAF4L2 WDR21C
Observed MW:	44 kDa
Uniprot ID :	Q8NA75

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09487(WDR21C antibody) at dilution of 1:500