

WDR6 antibody

Product Information

Catalog No.:	FNab09499
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Enhances the STK11/LKB1-induced cell growth suppression activity. Negative regulator of amino acid starvation-induced autophagy.

Immunogen information

Immunogen:	WD repeat domain 6
Synonyms:	tRNA (34-2'-O)-methyltransferase regulator WDR6 WD repeat-containing protein 6 WDR6
Observed MW:	122 kDa
Uniprot ID :	Q9NNW5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09499(WDR6 antibody) at dilution of 1:500