

WDR76 antibody

Product Information

Catalog No.:	FNab09502
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Specifically binds 5-hydroxymethylcytosine(5hmC), suggesting that it acts as a specific reader of 5hmC.

Immunogen information

Immunogen:	WD repeat domain 76
Synonyms:	WD repeat-containing protein 76 WDR76
Observed MW:	65 kDa
Uniprot ID :	Q9H967

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09502(WDR76 Antibody) at dilution of 1:300