

## XKR6 antibody

### Product Information

Catalog No.:	FNab09544
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	XK, Kell blood group complex subunit-related family, member 6
Synonyms:	XK-related protein 6 XKR6 C8orf21 C8orf5 C8orf7 XRG6
Observed MW:	42 kDa
Uniprot ID :	Q5GH73

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:5000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09544(XKR6 Antibody) at dilution of 1:1000