

## NHEJ1 antibody

### Product Information

Catalog No.:	FNab09545
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	nonhomologous end-joining factor 1
Synonyms:	Non-homologous end-joining factor 1 Protein cernunnos XRCC4-like factor NHEJ1 XLF
Observed MW:	33 kDa
Uniprot ID :	Q9H9Q4

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000
Image:	



human small intestine tissue were subjected to SDS PAGE followed by western blot with FNab09545(NHEJ1 antibody) at dilution of 1:300