

## YPEL3 antibody

### Product Information

Catalog No.:	FNab09569
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in proliferation and apoptosis in myeloid precursor cells.

### Immunogen information

Immunogen:	yippee-like 3(Drosophila)
Synonyms:	Protein yippee-like 3 YPEL3
Observed MW:	14 kDa
Uniprot ID :	P61236

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:5000
Image:	



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab09569(YPEL3 antibody) at dilution of 1:1000