

ZBTB44 antibody

Product Information

Catalog No.:	FNab09595
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger and BTB domain containing 44
Synonyms:	Zinc finger and BTB domain-containing protein 44 BTB/POZ domain-containing protein 15 Zinc finger protein 851 ZBTB44 BTBD15 ZNF851
Observed MW:	65 kDa
Uniprot ID :	Q8NCP5

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab09595(ZBTB44 Antibody) at dilution of 1:600