

ZG16 antibody

Product Information

Catalog No.:	FNab09632
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in protein trafficking. May act as a linker molecule between the submembranous matrix on the luminal side of zymogen granule membrane(ZGM) and aggregated secretory proteins during granule formation in the TGN.

Immunogen information


Immunogen:	zymogen granule protein 16 homolog(rat)
Synonyms:	Zymogen granule membrane protein 16 (Zymogen granule protein 16, hZG16) Secretory lectin ZG16 ZG16
Observed MW:	18 kDa
Uniprot ID :	O60844

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas using FNab09632(ZG16 antibody) at dilution of 1:100



human colon tissue were subjected to SDS PAGE followed by western blot with FNab09632(ZG16 antibody) at dilution of 1:500