

ZNF114 antibody

Product Information

Catalog No.:	FNab09654
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger protein 114
Synonyms:	Zinc finger protein 114 ZNF114
Observed MW:	45 kDa
Uniprot ID :	Q8NC26

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:200-1:2000
Image:	



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab09654(ZNF114 Antibody) at dilution of 1:600