

ZNF215 antibody

Product Information

Catalog No.:	FNab09668
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger protein 215
Synonyms:	Zinc finger protein 215 BWSCR2-associated zinc finger protein 2 (BAZ-2) Zinc finger protein with KRAB and SCAN domains 11 ZNF215 BAZ2 ZKSCAN11
Observed MW:	60 kDa
Uniprot ID :	Q9UL58

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:200-1:2000; IF: 1:10-1:100
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab09668(ZNF215 Antibody) at dilution of 1:300