

ZSCAN31 antibody

Product Information

Catalog No.:	FNab09685
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	zinc finger protein 323
Synonyms:	Zinc finger and SCAN domain-containing protein 31 Zinc finger protein 323 ZSCAN31 ZNF310P ZNF323
Observed MW:	47 kDa
Uniprot ID :	Q96LW9

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab09685(ZNF323 Antibody) at dilution of 1:600