

ZSCAN5A antibody

Product Information

Catalog No.:	FNab09759
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information

Immunogen:	zinc finger and SCAN domain containing 5A
Synonyms:	Zinc finger and SCAN domain-containing protein 5A Zinc finger protein 495 ZSCAN5A ZNF495 ZSCAN5
Observed MW:	56 kDa
Uniprot ID :	Q9BUG6

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:5000
Image:	



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab09759(ZSCAN5A Antibody) at dilution of 1:1000