
putrescine antibody

Product Information

Catalog No.:	FNab09929
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	putrescine
Synonyms:	1,4-Diaminobutane Putrescine dihydrochloride
Observed MW:	None
Uniprot ID :	None

Application

Reactivity:	None
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human intestinal cancer tissue slide using FNab09929 (putrescine antibody) at dilution of 1:200