

## TFAP2A antibody

### Product Information

Catalog No.:	FNab09980
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6C5
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	transcription factor AP-2 alpha
Synonyms:	Transcription factor AP-2-alpha (AP2-alpha) AP-2 transcription factor Activating enhancer-binding protein 2-alpha Activator protein 2 (AP-2) TFAP2A AP2TF TFAP2
Observed MW:	50 kDa
Uniprot ID :	P05549

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:2000 - 1:10000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab09980(TFAP2A antibody) at dilution of 1:5000