

anti- ACTN3 antibody

Product Information

Catalog No.: FNab10008

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Immunogen information

Immunogen: actinin, alpha 3

Synonyms: Actinin, alpha 3, ACTN3, Alpha actinin 3, F actin cross linking protein

Observed MW: 103 kDa Uniprot ID: Q08043

Application

Reactivity: Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:2000; IHC: 1:50-1:500

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded rat heart tissue slide using FNab10008(ACTN3 Antibody) at dilution of 1:200



rat skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab10008(ACTN3 Antibody) at dilution of 1:1000