

NABP2 antibody

Product Information

Catalog No.:	FNab10167
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	SOSS complex subunit B1
Synonyms:	SOSS complex subunit B1 Nucleic acid-binding protein 2 Oligonucleotide/oligosaccharide-binding fold-containing protein 2B Sensor of single-strand DNA complex subunit B1 Sensor of ssDNA subunit B1 (SOSS-B1) Single-stranded DNA-binding protein 1 (hSSB1) NABP2 OBFC2B SSB1 SOSB1
Observed MW:	22 kDa
Uniprot ID :	Q9BQ15

Application

Reactivity:	Human, rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded Rat testis tissue slide using FNab10167(SOSB1 antibody) at dilution of 1:200