

PAX4 antibody

Product Information

Catalog No.:	FNab10171
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Paired box protein Pax-4
Synonyms:	Paired box protein Pax-4 PAX4
Observed MW:	38 kDa
Uniprot ID :	O43316

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human kidney cancer tissue slide using FNab10171(PAX4 antibody) at dilution of 1:200