

MYBPC3 antibody

Product Information

Catalog No.:	FNab10193
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Myosin-binding protein C, cardiac-type
Synonyms:	Myosin-binding protein C, cardiac-type (Cardiac MyBP-C) C-protein, cardiac muscle isoform MYBPC3
Observed MW:	141 kDa
Uniprot ID :	Q14896

Application

Reactivity:	Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded mouse heart tissue slide using FNab10193(MYPC3 antibody) at dilution of 1:200