

SOS2 antibody

Product Information


Catalog No.:	FNab10316
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Son of Sevenless Homolog 2
Synonyms:	Son of sevenless homolog 2 (SOS-2) SOS2
Observed MW:	152 kDa
Uniprot ID :	Q07890

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab10316(SOS2 antibody) at dilution of 1:100



A549 cells were subjected to SDS PAGE followed by western blot with FNab10316(SOS2 Antibody) at dilution of 1:1000