

ITCH antibody

Product Information


Catalog No.:	FNab10329
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	E3 ubiquitin-protein ligase Itchy homolog
Synonyms:	E3 ubiquitin-protein ligase Itchy homolog (Itch) Atrophin-1-interacting protein 4 (AIP4) HECT-type E3 ubiquitin transferase Itchy homolog NFE2-associated polypeptide 1 (NAPP1) ITCH
Observed MW:	103 kDa
Uniprot ID :	Q96J02

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNab10329(ITCH antibody) at dilution of 1:100



THP-1 cells were subjected to SDS PAGE followed by western blot with FNab10329(ITCH antibody) at dilution of 1:1000