

FUT2 antibody

Product Information

Catalog No.:	FNab10367
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	Galactoside alpha-(1,2)-fucosyltransferase 2
Synonyms:	Galactoside alpha-(1,2)-fucosyltransferase 2 Alpha(1,2)FT 2 Fucosyltransferase 2 GDP-L-fucose:beta-D-galactoside 2-alpha-L-fucosyltransferase 2 SE2 Secretor blood group alpha-2-fucosyltransferase Secretor factor (Se) Type 1 galactoside alpha-(1,2)-fucosyltransferase FUT2 Type 2 galactoside alpha-(1,2)-fucosyltransferase FUT2 FUT2 SEC2
Observed MW:	40 kDa
Uniprot ID :	Q10981

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10367(FUT2 antibody) at dilution of 1:1000