

NAT2 antibody

Product Information


Catalog No.:	FNab10375
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Arylamine N-acetyltransferase 2
Synonyms:	Arylamine N-acetyltransferase 2 Arylamide acetylase 2 N-acetyltransferase type 2 (NAT-2) Polymorphic arylamine N-acetyltransferase (PNAT) NAT2 AAC2
Observed MW:	33 kDa
Uniprot ID :	P11245

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded rat lung slide using FNab10375(ARY2 Antibody) at dilution of 1:100



mouse kidney were subjected to SDS PAGE followed by western blot with FNab10375(ARY2 antibody) at dilution of 1:1000