

ITGAL antibody

Product Information


Catalog No.:	FNab10427
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	integrin, alpha L
Synonyms:	Integrin alpha-L CD11 antigen-like family member A Leukocyte adhesion glycoprotein LFA-1 alpha chain (LFA-1A) Leukocyte function-associated molecule 1 alpha chain ITGAL CD11A
Observed MW:	175 kDa
Uniprot ID :	P20701

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab10427(CD11A antibody) at dilution of 1:500