

BAZ1A antibody

Product Information


Catalog No.:	FNab10639
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Bromodomain adjacent to zinc finger domain protein 1A
Synonyms:	ACF1 BAZ1A CHRAC subunit ACF1 DKFZp586E0518 FLJ14383 hACF1 hWALp1 WALp1 WCRF180
Observed MW:	250 kDa
Uniprot ID :	Q9NRL2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10639(BAZ1A Antibody) at dilution of 1:1000