
GLI1 antibody

Product Information


Catalog No.:	FNab10657
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	Zinc finger protein GLI1
Synonyms:	GLI, GLI family zinc finger 1, GLI1, Glioma associated oncogene, Oncogene GLI, Zinc finger protein GLI1
Observed MW:	150 kDa
Uniprot ID :	P08151

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded rat lung using FNab10657(GLI1 antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab10657(GLI1 antibody) at dilution of 1:1000