
PSA antibody

Product Information

Catalog No.:	FNab10670
Size:	1mg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1G9
IsoType:	IgG1
Storage:	PBS with 0.1% Proclin300 pH 7.3, store at 2-8° C up to 21 days, and at -20° C for long-term storage. (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	kallikrein-related peptidase 3
Synonyms:	APS, Gamma seminoprotein, hK3, Kallikrein 3, kallikrein related peptidase 3, KLK2A1, KLK3, KLK3, PSA, P 30 antigen, Prostate specific antigen, PSA, Semenogelase, Semin
Uniprot ID :	P07288

Application

Reactivity:	Human
Tested Application:	ELISA, CLIA (detection antibody)
Recommended dilution:	Each laboratory should determine an optimum working titer for use in its particular application.