
OSBPL3 antibody

Product Information


Catalog No.:	FNab10671
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	OSBP related protein 3
Synonyms:	KIAA0704, ORP 3, ORP3, OSBP related protein 3, OSBP3, OSBPL3
Observed MW:	100 kDa
Uniprot ID :	Q9H4L5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab010671(OSBPL3 antibody) at dilution of 1:100



HeLa lysates were subjected to SDS PAGE followed by western blot with FNab10671(OSBPL3 antibody) at dilution of 1:1000