
SPC24 antibody

Product Information


Catalog No.:	FNab10681
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Kinetochose protein Spc24
Synonyms:	hSpc24, Kinetochose protein Spc24, SPBC24, SPC24
Observed MW:	26 kDa
Uniprot ID :	Q8NBT2

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Hela cells were subjected to SDS PAGE followed by western blot with FNab10681(SPC24 antibody) at dilution of 1:1000