

## HAS1 antibody

### Product Information


Catalog No.:	FNab10688
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Hyaluronan synthase 1
Synonyms:	Hyaluronan synthase 1 HA synthase 1  HuHAS1  HAS1 HYAS1
Observed MW:	65 kDa
Uniprot ID :	Q92839

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



Various lysates were subjected to SDS PAGE followed by western blot with FNab10688(HAS1 antibody) at dilution of 1:1000