

## CDH6 antibody

### Product Information

Catalog No.:	FNab10805
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background


This gene encodes a member of the cadherin superfamily. Cadherins are membrane glycoproteins that mediate homophilic cell-cell adhesion and play critical roles in cell differentiation and morphogenesis. The encoded protein is a type II cadherin and may play a role in kidney development as well as endometrium and placenta formation. Decreased expression of this gene may be associated with tumor growth and metastasis.

### Immunogen information

Immunogen:	K-Cadherin
Synonyms:	Cadherin-6 Kidney cadherin (K-cadherin) CDH6
Observed MW:	100 kDa
Uniprot ID :	P55285

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500 - 1:2000
Image:	



mouse brain were subjected to SDS PAGE  
followed by western blot with FNab10805(CDH6  
antibody) at dilution of 1:1000